# SILAVEST GOLD

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## DIN EN ISO 15912

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### DIN EN ISO 15912

#### Description:

Silavest Gold is a phosphate-bonded, graphite-free speed cast investment with exceptional properties. It is suitable for casting all types of precious alloys for the crown and bridge work technique.

#### Waxing up:

Optimum fit and surface quality are best achieved using patterns fabricated entirely from wax. Using resins often results in much poorer surface quality and cracks in the mould with the speed cast technique. This applies particularly to waxing up removable sections of attachment and telescopic crown work.

#### Technical data:

|                                  |   |
|----------------------------------|---|
| Mixing liquid:                   | SILADENT type 100 expansion liquid                |
| Mixing ratio<br>Powder : liquid: | 100 g : 22 ml                                     |
| Mixing under vacuum:             | 60 sec.   |
| Working time:                    | 5 - 7 Min.  |
| Storage and working temperature: | 17 - 19 °C  |
| Initial setting time:            | 9 - 11 min.                                       |
| Setting expansion:               | approx. 1.00%<br>(100% undiluted type 100 liquid) |

#### Attaching sprues:

Follow the instructions in the SILADENT-system when sprueing wax wires especially with Perawax sprues. The patterns should never be placed in the heat centre of the mould. Thin sections, e.g. crown margins, should face towards the sides of the mould and be 5 mm from the casting ring. SILADENT silicone sleeves produce a special surface structure of the moulds, which gives the investment a constant stiffness. For anterior tooth bridges we recommend to use our round silicone sleeves size 6 (REF 102610), for back tooth bridges resp. circularly constructions our model form silicone sleeves sizes 0 - 3.

#### Concentration of the expansion liquid:

22 ml liquid / 100 g powder, 35 ml liquid / 160 g sachet

| Concentration (%) | 1 x 160 g sachet |             | 2 x 160 g sachet |             | 3 x 160 g sachet |             |
|-------------------|------------------|-------------|------------------|-------------|------------------|-------------|
|                   | liquid (ml)      | dist. water | liquid (ml)      | dist. water | liquid (ml)      | dist. water |
| 100               | 35               | -           | 70               | -           | 105              | -           |
| 95                | 33               | 2           | 66               | 4           | 99               | 6           |
| 90                | 32               | 3           | 64               | 6           | 96               | 9           |
| 85                | 30               | 5           | 60               | 10          | 90               | 15          |
| 80                | 28               | 7           | 56               | 14          | 84               | 21          |
| 75                | 26               | 9           | 52               | 18          | 78               | 27          |
| 70                | 25               | 10          | 50               | 20          | 75               | 30          |
| 65                | 23               | 12          | 46               | 24          | 69               | 36          |
| 60                | 21               | 14          | 42               | 28          | 63               | 42          |
| 55                | 19               | 16          | 38               | 32          | 57               | 48          |
| 50                | 17.5             | 17.5        | 35               | 35          | 52.5             | 52.5        |
| 45                | 16               | 19          | 32               | 38          | 47               | 58          |

### Expansion control %:

|                                      | Casting alloys<br>high gold<br>content | Casting alloys<br>reduced<br>gold content | Bonding alloys<br>high gold<br>content | Bonding alloys<br>reduced<br>gold content |
|--------------------------------------|--|---|--|---|
| Inlays, onlays, partial crowns (wax) | 50 - 75 ± 5                            | 50 - 75 ± 5                               | -                                      | -   |
| Crown and - bridgework (wax)         | 65 ± 5                                 | 65 ± 5                                    | 85 ± 5                                 | 85 ± 5                                    |
| Telescope crowns (pattern resin)     | 75 ± 5                                 | 75 ± 5                                    | 80 ± 5                                 | 80 ± 5                                    |

### Mixing:

Use a spatula to mix the investment material by hand until the powder has been wetted thoroughly. Mix the investment for 60 seconds. Check the vacuum unit frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the surface. The working time from the investment is approx. 6 minutes.

### Investing:

Fill the ring with investment. Only use a vibrator if it is required to improve the flow. Avoid high vibration! This will lead to the formation of bubbles and breakdown of the mixture.

### Speed preheating:

Place the ring after 30 minutes in the max. 850°C preheated furnace. If a preheating temperature higher than 850°C is required, the temperature can be adjusted 20 minutes after placing the mould in the furnace.

**Caution:** During the first 15 minutes the furnace door should not be opened because of the risk of instant combustion.

### Programmed preheating:

|                 | Temperature          | Heating rate | Holding times |
|-----------------|----------------------|--------------|---------------|
| 1. Heating rate | Room temp. - 270 °C  | 8 °C / min.  | 30-50 min.    |
| 2. Heating rate | 270 °C - 580 °C      | 8 °C / min.  | 30 min.       |
| 3. Heating rate | 580 °C - Final temp. | 8 °C / min.  | 40 min.       |

### Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

### Caution:

**Inhaling silica dust is a health risk. A suitable face mask should be worn.**

**Further questions:** Contact SILADENT technical services  
(Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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