# SILADENT

Dr. Böhme & Schöps GmbH

Das System

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Dr. Böhme & Schöps GmbH

The System

Model creation  
Duplicating  
Investing  
Casting

## TEK-1 LEG

DIN EN ISO 22674

€ 0124

Composition: (in % by mass)	Technical data: (guidelines)
Co 58.0	Proof stress $R_{p0,2}$ (MPa) 480
Cr 27.0	Elongation at rupture $A_5$ (%) 10
Mo 6.0	Modulus of elasticity (GPa) 200
W 6.0	Coefficient of expansion
	25°C-500°C $14.0 \times 10^{-6}$
	25°C-600°C $14.2 \times 10^{-6}$
<b>Other constituents:</b> Si, Fe, Mn	Tensile strength $R_m$ (MPa) 480
Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.	Vickers hardness VH 10 285
	Density (g/cm <sup>3</sup> ) 8.4
	Melting range (°C) 1.270-1.380

### Description:

TEK-1 LEG is a non-precious, chromium cobalt-based bonding alloy. It does not contain any nickel, beryllium or gallium. TEK-1 LEG provides an economical alternative to precious metal and palladium-based alloys.

Its low hardness of only 285 VH 10 allows TEK-1 LEG to be easily milled, which is an advantage when fabricating fixed units for fixed/removable appliances.

### Investing and casting:

TEK-1 LEG can be used with phosphate-bonded investments, such as TEK-1 Vest, Presto Vest II, Premium and TeleVest. It is preheated to 850 °C - 950 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the SILADENT Technique instructions when attaching sprues. Use a separate ceramic crucible for casting TEK-1 LEG high-grade, non-precious alloy to prevent it being contaminated by other alloys. Clean the crucible after each cast. Cylindrical, continuous cast ingots collapse more easily during induction melting of the alloy. The alloy is cast when all the cylindrical ingots have melted together and a **final shadow runs over the molten metal just before the oxide layer disintegrates**. Please note that the melting point may vary depending on the type of casting machine used and when melting with vacuum, for example, less oxide forms and the oxide layer in the vacuum disintegrates more quickly. Considerable experience is required when melting high-grade alloys automatically, and it is most important that heating is set exactly according to the manufacturer's instructions. Inaccurate heat setting impairs the alloy and can have a detrimental effect on the porcelain bond. Always use the reducing zone of the flame in a circular motion to melt the alloy. Casting TEK-1 LEG repeatedly is not recommended, as this could alter both its composition and properties.

### Firing porcelain:

Standard porcelains, such as those manufactured by Vita, Ivoclar, DeTrey, Dentsply or Ducera, can be used.

Follow the relevant manufacturer's instructions for use.

1. **Sandblast** the surface in a pencil blaster with **150 µm - 250 µm**.
2. **Clean** the framework with distilled water or bioalcohol in an ultrasonic cleaner.
3. Carry out **oxide firing for 5 minutes at 960 °C - 980 °C** with vacuum.
4. **Always sandblast the oxide layer again after the oxide firing and clean it again.**
5. Apply the **first opaquer firing** thinly. Only cover uniformly at the second opaquer firing.
6. Fire according to the porcelain manufacturer's instructions. Clean the framework before each firing.
7. We recommend **slow cooling** for each dentine, second and glaze firing with a cooling phase to approx. 600°C.

### Soldering and welding:

Presolder using SILADENT Cr-Co solder and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys. Laser weld using a SILADENT Cr-Co laser welding rod (REF 102806).

TEK-1 LEG is classified as a class IIa medical product according to the European directive 93/42/EWG. SILADENT Dr. Böhme & Schöps GmbH packaging carries the € 0124 mark.

**Further questions:** Contact SILADENT technical services  
(Tel.: +49 (0) 53 21 - 37 79 25 / 26)  
or our sales representatives.

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