TEK-1 Vest

DIN EN ISO 15912

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DIN EN ISO 15912

Description:

TEK-1 Vest is a phosphate-bonded, graphite-free precision investment for TEK-1 one-piece-casting and for partial denture framework.

Casting models:

Measure the required amount of investment and the correct concentration of expansion liquid according to the mixing ratio for casting models. Mix manually for 20 sec. with a spatula before vacuum mixing for 60 sec. Use small vibration for pouring the contours. An instrument is useful for this. Switch off the vibrator and pour the rest of the model.

Use of pressure:

Investing under pressure is **not necessary** due to TEK-1 Vest's extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 8 - 10 minutes so that it does not impede the setting expansion.

Concentration of the expansion liquid: 20 ml Liquid / 100 g powder

Mixing ratio:

(Depending on the laboratory, the specifications may vary individually.) The examples apply for silicone duplications. Mixing ratio powder/liquid: 100 g powder : 20 ml liquid.

By the **TEK-1 one-piece-casting** we recommend a two-stage investment: For telescopes, depending on the friction 60 % - 70 %. In the ideal case 65 % (13 ml liquid and 7 ml distilled water per 100 g powder). After that the model gets effused with a liquid concentration of 55 % by using the wet method.

50 g powder	6,5 ml liquid	3,5 ml distilled water	For the crowns 65 %
150 g powder	16,5 ml liquid	13,5 ml distilled water	For the modell 55 %
400 g powder	44,0 ml liquid	36,0 ml distilled water	For the mould 55 %

Framework with clasps	150 g Pulver	16,5 ml	13,5 ml	55 %
Fixed / removable work	150 g Pulver	16,5 ml	13,5 ml	55 %

By **fixed-removable work** we recommend 65 % - 75 % liquid concentration according to the size and the expansion of the jaw. By **Galvano tertiar framework** we recommend to manufacture the

Technical data:

Mixing liquid:	SILADENT TEK-1 liquid
Mixing ratio, Powder : liquid:	100 g : 20 ml
Mixing under vacuum:	60 sec.
Working time:	4 - 5 min.
Recommended storage/ Working temperature:	17 - 19°C
Initial setting time:	7 - 10 min.
Preheating temperature range for placing in the furnace:	900 - 1.000 °C
Place in the hot furnace:	25 min. after investing

two-stage investment by using 100 % liquid concentration for the frame and 55 % for the model.

Removing the duplicate model:

Carefully remove the investment model after approx. 20 min. Start waxing up once the model is cool enough to handle and the wax is no longer liable to melt.

Investing:

Measure and mix TEK-1 Vest as described above. For the investing use the same liquid concentration as for the model. A small reduction can optimize the adhesion between model and investing. Pour with minimum vibration until the pattern is completely covered. Switch off the vibrator, fill the ring and do not vibrate again. Moulds intended for speed casting should be placed on a very damp cloth so that the dehydrated model regains some moisture and achieves a better bond with the investment in the mould.

Speed casting:

Exact 25 minutes after investing, place the mould with the sprue hole pointing downwards in a furnace preheated to the correct temperature (900 - 1.000 °C). We recommend a maximum preheating temperature of 900 °C. This temperature is perfectly adequate for both casting the mould and final thermal expansion. Moulds are ready to cast approx. 60 min. after the final temperature is reached. During the first 30 min. after inserting the moulds, the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion.

Conventional preheating:

Place the mould into the cold furnace after 25 minutes.

	Temperature	Heat rate	Holding stages
1. heat rate	280 °C	3 - 7 °C / minute	45 - 60 minutes
2. heat rate	580 °C	3 - 7 °C / minute	45 - 60 minutes
3. end temperature	900 °C	3 - 7 °C / minute	30 - 45 minutes

Maximum end temperature: 1.050 °C when investing overnight the expansion is by 5 % larger. This should be considered when deciding about the liquid concentration.

Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

Caution:

Inhaling silica dust is a health risk. A suitable face mask should be worn.

Further questions: Contact SILADENT technical services (Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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Dr. Böhme & Schöps GmbH